

# DC-8

Diagnostic Ultrasound System



confidence | comfort | ease



## Introducing the DC-8: Designed by users for users

Healthcare continues to evolve and new tools are developed to improve clinical functions and simplify processes. Our products are continually refined to ease the tasks of clinicians and increase patient throughput. The DC-8, designed by users for users, was developed in the spirit of this commitment.





The DC-8 Ultrasound System supports the most demanding applications on a day-to-day basis. Advanced technologies, superior image quality, and a special focus on ergonomic design combine to improve patient management, increase ease of use, and enhance users' overall level of confidence.



Intuitive touch screen



Adjustable control panel



Articulating arm adjusts height and angle

DC-8's **ergonomic solutions** address every part of your body that interacts with the device. We know that users have a physical relationship with the device and with this fact in mind, we have improved procedural efficiency and reduced repetitive actions. A new user interface enables logical and efficient workflow while ensuring more convenient operation.

The result: easier acquisition of clinical images acquired with less stress. What's more, the intuitive touch screen provides direct access to advanced modes and features. Working together, these features and benefits simplify and enhance the clinical experience.

#### **User Centric Control Panel**

- Common commands are easily executed with minimal movement of one hand
- The height-adjustable control panel swivels and rotates effortlessly
- Control keys, dials, and toggles, intuitively grouped within an innovative kidney-shape design, enhance comfort and support healthy posture

#### **Ergonomic Solutions**

- Lightweight transducers designed with soft flexible cables reduce stress while scanning
- 10.4" touch panel places all functions within reach
- 19" high-definition LCD monitor helps reduce eye fatigue and can be viewed at angles greater than 170°
- Wide degree of monitor tilt and rotation works with your natural body movement to ensure maximum comfort
- Monitor fold-down function allows for easy mobility

## See more

DC-8's large 19" high-definition LCD monitor provides better image detail.

## Save time

iWorks/iWorks OB reduces the number of key strokes required for a complete diagnostic study by 71-82%.

## Eliminate tasks

The user interface reduces repetitive motions required to complete complex ultrasound functions.

## Move gracefully

Height-adjustable control panel swivels effortlessly.

## Reduce stress

The system design exceeds industry standards to reduce work related injury.






When we see,



Sclerosed CCA | L14-6WE



Femoral artery bifurcation | L14-6WE




Liver with ascites and pleural effusion | C5-2E



Kidney perfusion | C5-2E




Adult heart displaying MR and TR | P4-2E




TDI with quantitative analysis | P4-2E


we understand.



Thyroid mass | L14-6WE



Testicular cyst | L14-6WE



Fetal face (30 week) surface rendering mode | 4D probe



Fetal profile (14-week) surface rendering mode | 4D probe



Ovarian mass | V11-3E



Fetal cord | C7-3E



Fetus profile using 3T technology



Thyroid mass using HR Flow

## Clarity = better outcomes

**Multi-Beam Parallel Imaging Technology™** integrates advanced imaging technologies to produce optimized images with minimal user involvement. When combined with the DC-8's echo-enriched beamformer, it produces excellent images.

**New 3T Transducers** use Mindray's patent-pending design, which increases bandwidth and transmission efficiency, providing better resolution and S/N ratio.

- Scans with higher frequencies even at deeper depths
- Captures images easier and faster in different body organs
- Increased image contrast and spatial resolution while scanning the most challenging patients

**HR Flow** is an innovative technology that is designed to visualize vessel profiles and complex flow patterns.

- Provides detection of color Doppler signals
- Offers extraordinary spatial resolution without compromising penetration or sensitivity





### Advanced Imaging Technologies

- **iBeam™** spatial compounding using multiple scan angles to form a single image with enhanced contrast resolution for better visualization
- **iScape™** provides real-time panoramic imaging of large structures in a single sweep
- **iTouch™** offers single-button image optimization that provides excellent image uniformity
- **iClear™** speckle suppression and edge enhancement technology improves contrast resolution to differentiate between normal and abnormal structures
- **iZoom™** one-touch expands displayed image for viewing at a distance
- **iNeedle™** enhances needle visualization to aid in biopsy procedures





DC-8

IWorks

Exit

Repeat	Replace	Auto Calc.	Previous	Next
iNCert	Delete	Invert	Stop	Suspend

Quality	Persistence	Dyn Range	Gray Map	Tint
Harmonic On				On
General	7 / 8	240 dB	6 / 8	6 / 8

A. Power

Volume



## More intelligent workflow

DC-8 incorporates tools to increase your speed, accuracy, and productivity while maintaining the flexibility you need to complete exams comfortably and efficiently.

**iWorks** is a powerful protocol management system beneficial for all ultrasound practitioners. From an experienced sonographer to a student, or the occasional visiting physician, the DC-8 allows the user to focus on the patient and spend less time on operational steps and keystrokes. The system:

- Automatically guides the user through an exam protocol
- Automatically inserts comments
- Automatically sends measurement results to the report for easy interpretation
- Enables customized protocols for unusual cases
- iNSert is a unique feature that allows users to measure and document common pathology within iWorks

**iWorksOB** is a dedicated extension of iWorks, designed to work easily with the random movements of the fetus.

- Practical solution for OB automated protocols
- Automatically sends measurement results to the report for easy interpretation
- Recognizes the random nature of OB exams to better facilitate workflow

**Raw Data Manipulation** is a versatile and powerful tool for optimizing images and analyzing information even after the image has been saved. By providing this flexibility, the DC-8 enables the user to work with the stored data wherever and whenever it is needed.

- Allows you to manipulate multiple parameters, measurements, and annotate on acquired images

**Seamless connectivity.** Plug and play offers you a wireless LAN option to connect to your local network wherever you are. Transfer images or patient data with ease and convenience. Comprehensive DICOM features: Store/Print, Worklist, MPPS, Query/Retrieve, OB/GYN S/R, Vascular S/R and Cardiac S/R.



Mindray global offices

## Proceed with confidence

Mindray is a global medical device business with more than 6,000 employees, built on three established product segments: Patient Monitoring and Life Support Products, *In-Vitro* Diagnostic Products, and Medical Imaging Systems.

We are proud to manufacture 95% of the components used in our products. Our in-house capabilities ensure quality control, rapid time-to-market, and reduced manufacturing and transportation costs. These competencies enable Mindray to provide clinicians with highly regarded quality solutions at an outstanding value. These products and services are at the heart of the commitment we make to both our customers and, ultimately, the patients in their care.

In addition, Mindray's clinical education specialists provide your staff with comprehensive, on-site training and service. Our dedicated field service and in-house technical support organizations are among the largest direct-service teams in the industry. Individually and collectively they are committed to enhancing patient care while protecting your investment. Learn how we can help you maximize your resources with some of the most advanced and affordable devices in the industry. For more information about Mindray products, visit [www.na.mindray.com](http://www.na.mindray.com).



Mindray North America  
U.S. headquarters in Mahwah, NJ

Mindray DS USA, Inc.  
800 MacArthur Blvd., Mahwah, NJ USA 07430  
Tel: 1.800.288.2121 or 201.995.8000 Fax: 1.800.926.4275 [www.na.mindray.com](http://www.na.mindray.com)

